

Modelos "Housing Fish"

Xenopus and fish systems

- Extra large window for a clear view of the fish.
- Equipped with four filtration systems to maintain a very good quality of the water.
- Low water alarm for the water tank.
- Shelving made of AISI 316 stainless steel, RESISTANT TO SALT WATER EROSION.
- Tanks made of PVC, RESISTANT TO ACIDS AND BASES. .
- Independent water regulation valve.
- UV sterilization ensures bacteria-free water.



Ordering Information

Reference	Description	Dimensions
67-AD-026-A	SIST. AUTONOMOUS 2 FLOORS FOR XENOPUS AND FISH; INCLUDES: AIR/WATER PUMPS, PRE-FILTER, FILTERS: BAG+ CARBON+ BIOLOGICAL. UV STERILIZATION	142 x 47 x 190 cm